


# X-Dock™

Docking station for X-Series™ instruments

▶ The new X-Series™ instruments are Alpha Technologies' modular production floor solution for customers who are looking to bring testing closer to the production line.



▶ The X-Series™ brings Alpha's superior testing solutions to the production floor. The "X" in the X-Series™ brand name stands for "cross-over" from lab environment to production floor.



**ALPHA**TECHNOLOGIES

## Features




- Accommodates up to 14,000 tests between film roll changes
- Touch screen industrial computer
- Configurable to use 2, 3, or 4 layers of film
- Organize data into reports and perform statistical analysis with integrated Enterprise software
- Multiple languages available



## Performance



- Compatible with MDR-X™ and RPA-X™
- Auto test module alignment
- Data management solution
- Production floor certified safety compliance
- Network ready
- Test configuration management



## Specifications



DIMENSIONS: W: 27.6 in (70 cm), D: 28.9 in (73.4 cm), H: 61.3 in (155.7 cm)      ELECTRICAL: 85-250 VAC/20A, 50-60 Hz

