

densiTECH+™

Automatic density measuring instrument designed for rapid and accurate measurements of density. An American-made product available worldwide with good performance at a reasonable cost.

Features

- Analytical balance (available for other measurements).
- Simplified operation.

Specifications

- Tabletop unit.
- Balance measures to nearest 0.0001 g.

Performance

- Measures the density of up to four samples per minute.
- Accuracy is within $\pm 0.025\%$.
- Balance resolution is good to 0.0001 g.

Test Capabilities

- Measure the density of cured rubber compounds.



densiTECH+™

www.alpha-technologies.com